

## Infrared Detector

### USE OF COMBINED GC/IRD/MSD IN THE CHARACTERIZATION OF SYNTHETIC BUTTER FLAVOR

#### Introduction

There are several types of synthetic butter flavors used today. Some contain mostly diacetyl with vanillan, various butyrates, and butyric acid. There are other types of synthetic butter flavors also. One of those was analyzed here. It surprisingly contains non of the expected smaller flavorants but higher molecular weight lactones.

#### Equipment

For the combined GC/IRD/MSD system, the gas chromatograph, an HP 5890 Series II was set up using a 30 meter DB-5 column. The column effluent was split at the end of the column at a 10 to 1 ratio with the bulk of the flow going to the HP 5965B IRD and the lesser amount going to the HP 5971A MSD. The details of this parallel configuration have been previously described.

Control of the GC, IRD and MSD as well as data reduction was done with G1094A software on the single PC based G1030A MS/IR ChemStation (DOS).

#### Results

With the combined GC/IRD/MSD high confidence qualitative analysis is possible. Using this sample of synthetic buter flavor several features of the system have been utilized, viz., library searching of infrared and mass spectra with combined reporting and selected wavenumber (functional group) chromatography.

Figure 1 is a portion of the IRD total response chromatogram (TRC). The peak at 12.18 minutes turns out to be delta dodecalactone, see Figure 2. The library search results for this peak are in Figure 3. IRD library searching allows various libraries to be chained and over 11,000 spectra were searched. Note in Figure 2 that names and structures can be annotated.

The combined IR and MS display is in Figure 4. Here the total response chromatogram and the total ion chromatogram have been software aligned for better presentation. The IR and Mass Spectra for the peak at about 12.18 minutes are also displayed.

ASAP, LLC  
1511 Neave Street  
Covington, KY 41011  
ph: 877-987-2800 or  
859-581-6990  
fx: 859-581-6821  
info@asapanalytical.com  
www.asapanalytical.com  
date founded: 1999





## Infrared Detector




Figure 4. Combined IRD and MSD chromatograms and spectra for peak at 12.18 min




Figure 6. IRD Total Response Chromatogram and Selected Wavenumber Chromatograms for Lactones

File: C:\PROGRAM FILES\BUTTERFLY.D  
 Operator: Roger Leland  
 Acquired: 12 Jun 02 3:17 pm Using AcqMethod BUTTER.M  
 Sample Name: Synthetic Buter Peak  
 Misc Info: 30M 250M SUI DBS 526 5920+750278201 82  
 Search Libraries: C:\DATABASES\NIST1 Minimum Quality: 95  
 C:\DATABASES\NIST75K Minimum Quality: 0

Unknown Spectrum: Acet  
 Integration Param: AutoIntegrate

RET	AREA%	LIBRARYID	REF	CASE	QUAL
1	9.70	C:\DATABASES\NIST1			
		2(3H)-Furanone, 5-hexylhydro	121907	000705-14-3	85
		2(3H)-Furanone, 5-hexylhydro	121908	000705-14-3	85
		GAMMA-BUTYROLACTONE, 5-HEXYL-	22695	000009-09-3	83
IR2	9.71	C:\DATABASES\REFL\REFL.LIB			
		SUCCINIC ANHYDRIDE, 11-OCTYL-11A	1921	007757-95-2	9648
		C:\DATABASES\REFL\REFL.LIB			
		GAMMA-OCTANOLACTONE, 5-HEXYL-	2241	000194-00-7	9963
		C:\DATABASES\REFL\REFL.LIB			
		UNDECANOL, GAMMA-LACTONE, 9-HEXYL-	2244	000104-07-0	9794
2	10.70	C:\DATABASES\NIST75K			
		2H-Pyran-2-one, 6-(4-hydroxy-5-pentyl)	88208	000705-06-2	60
		2H-Pyran-2-one, 6-(4-hydroxy-5-pentyl)	88223	000713-06-1	95
		1,3,2-Dioxaborolane, 3,4-dithienyl	4896	057698-69-3	93
IR3	10.73	C:\DATABASES\REFL\REFL.LIB			
		(E)-LACTIDIC ACID, DELTA-LACTONE, 9B	2259	000705-06-2	1792
		C:\DATABASES\REFL\REFL.LIB			
		(E)-LACTIDIC ACID, DELTA-LACTONE, 9B	2260	000710-04-3	1792
		C:\DATABASES\REFL\REFL.LIB			
		2H-PYRAN-2-ONE, 5-HEXYL-TETRAHYDRO-2H	1588	000710-04-3	2174

Figure 5. Combined IRD and MSD library search results for peak at 12.18 minutes

ASAP, LLC  
 1511 Neave Street  
 Covington, KY 41011  
 ph: 877-987-2800 or  
 859-581-6990  
 fx: 859-581-6821  
 info@asapanalytical.com  
 www.asapanalytical.com  
 date founded: 1999



## Infrared Detector

### Conclusion

As seen with this analysis of a synthetic butter flavor, the combined GC/IRD/MSD provides a powerful tool for the characterization of flavor components using library searching and selected wavenumber chromatography

### Conditions

#### Gas Chromatograph

*Column:* 30m x 0.25mm x .5um DB-5

*Carrier Gas:* Helium at 180kPa, 1mL/min

*Oven:* 50°C (2 min) at 250°C at 20°C/min with 2 min hold

*Injection Port:* 250°C

*Sample Injection:* 1ul, split 20 to 1

#### IRD Parameters

*Light Pipe:* 260°C

*Transfer Lines:* 280°C

*Sweep Gas:* Nitrogen 110 kPa inlet, 100 kPa outlet

*Scan Parameters:* 8cm<sup>-1</sup> resolution, 2 coads, 3 scans/second

*Detector:* Wide Band MCT, 4000-550cm<sup>-1</sup>

#### MSD Parameters

*Mass Range:* 10 to 310 daltons

*Scan Parameters:* 2 a/d samples, 1.4 scans/second stored

ASAP, LLC  
1511 Neave Street  
Covington, KY 41011  
ph: 877-987-2800 or  
859-581-6990  
fx: 859-581-6821  
info@asapanalytical.com  
www.asapanalytical.com  
date founded: 1999

