

Infrared Detector

WHAT'S NEW ON THE IRD³

The IRD 3 is the next generation in vapor phase infrared detectors from ASAP Analytical. The IRD 3 in combination with GC provides both separating and molecular identification in one system. The IRD 3 has new electronics, operating software and data analysis software that provides the user with increased performance and reliability.

SOFTWARE

- ASAP Stand Alone Instrument Control (ASIC)
 - Replaces IRD Overlay and G1701BA MSD ChemStation
 - Platform Independent
 - Windows XP and Windows 7 Compatible
- Essential FTIR (eFTIR)
 - Replaces GRAMS/32
 - Report Generation & Library Searching
 - View Spectra in 3-D
 - Selected Wavenumber Chromatogram (SWC)

HARDWARE

- Electronics
 - Eliminates Smart Card External Controller
 - Replaces GPIB with USB 2.0
 - Power Supply Board
 - Heater Control Board
 - Main CPU Board
 - Detector Preamp.
 - Laser Detector Board

ONE TECHNIQUE TO CONFIRM THEM ALL!



ASAP, LLC
1511 Neave Street
Covington, KY 41011
ph: 877-987-2800 or
859-581-6990
fx: 859-581-6821
info@asapanalytical.com
www.asapanalytical.com
date founded: 1999